

# Technical Specifications for the Cadmium Ion-Selective Electrode ELIT 8241

## Introduction

The Cadmium Ion-Selective Electrode has a solid-state crystal membrane. The electrode is designed for the detection of cadmium ions ( $\text{Cd}^{+2}$ ) in aqueous solutions and is suitable for use in both field and laboratory applications.

The Cadmium Ion is a divalent cation .

One mole of ( $\text{Cd}^{+2}$ ) has a mass of 112.411 grams; 1000ppm is 0.009M

Dissolve 2.744g Cadmium Nitrate tetra hydrate ( $\text{Cd}(\text{NO}_3)_2 \cdot 4\text{H}_2\text{O}$ ) in 1 Litre water.

## Physical Specification

Length of body excl gold contacts	130 mm
Length of body incl. gold contacts	140 mm
Diameter of body	8 mm
DC resistance at 25 °C	< 2.5 MOhm
Minimum feasible sample volume	5mls

## Chemical / Operational Specifications

Preconditioning / Standard solution <i>(But see General Operating Instructions)</i>	Normally 1000ppm $\text{Cd}^{+2}$ as $\text{Cd}(\text{NO}_3)_2$
Preconditioning time	5 minutes
Optimal pH range	pH 3 to pH 7
Temperature range	0 to 80 ° C
Recommended ISAB	5M $\text{NaNO}_3$ (add 2% v/v)
Recommended reference electrode	double junction ( <b>ELIT 003</b> )
Reference electrode outer filling solution	$\text{CH}_3\text{COOLi}$
Electrode slope at 25 °C	26±3 mV/ decade
Concentration range	0.1 to 11,000 ppm ( $9 \times 10^{-7}$ to 0.1 Molar)
Response time	< 10 seconds
<i>(Defined as time to complete 90% of the change in potential after immersion in the new solution.)</i>	
Potential drift <i>(in 1000 ppm)</i>	< 3 mV/ day (8 hours)
<i>(Measured at constant temperature and with ISE and Reference Electrode continually immersed)</i>	

## Interference:

NB: All poly-crystalline membranes contain Silver Sulphide and thus will not give reliable readings if Ag or S ions are present in the solution. Furthermore, this electrode cannot be used in the presence of significant concentrations of Copper or Mercury, and any Iron or Lead ions in a concentration greater than one hundredth of that of the Cadmium will lead to spuriously high values – Selectivity coefficients for Iron and Lead are approximately 10.  
Note low pH range for this electrode.

**For more information, see: [www.nico2000.net](http://www.nico2000.net).**