

# Technical Specifications for the Iodide Ion-Selective Electrode ELIT 8281

## Introduction

The Iodide Ion-Selective Electrode has a solid-state Crystal membrane. The electrode is designed for the detection of iodide ions (I<sup>-</sup>) in aqueous solutions and is suitable for use in both field and laboratory applications.

The Iodide Ion is a monovalent anion .

One mole of ( I<sup>-</sup>) has a mass of 126.905 grams; 1000ppm is 0.008M

Dissolve 1.308g anhydrous potassium iodide (KI) in 1 litre water.

## Physical Specifications

|                                    |           |
|------------------------------------|-----------|
| Length of body excl gold contacts  | 130 mm    |
| Length of body incl. gold contacts | 140 mm    |
| Diameter of body                   | 8 mm      |
| DC resistance at 25°C              | <2.5M Ohm |
| Minimum feasible sample volume     | 5mls      |

## Chemical / Operational Specifications

|  |  |
|--|--|
| Preconditioning/ Standard solution<br><i>(But see General Operating Instructions)</i>                    | Normally 1000ppm I <sup>-</sup> as KI                |
| Preconditioning time   | 5 minutes  |
| Optimal pH range   | pH 2 to pH 12  |
| Temperature range  | 0 to 80°C  |
| Recommended ISAB   | 5M NaNO <sub>3</sub> (add 2%v/v)                     |
| Recommended reference electrode  | double junction ( <b>ELIT 003</b> )                  |
| Reference electrode outer filling solution   | 0.1M CH <sub>3</sub> COOLi                           |
| Electrode slope at 25°C  | 54±5 mV/ decade                                      |
| Concentration range  | 0.06 to 12,700 ppm (5x10 <sup>-7</sup> to 0.1 Molar) |
| Response time  | < 10 seconds   |
| <i>(Defined as time to complete 90% of the change in potential after immersion in the new solution.)</i> |  |
| Potential drift <i>(in 1000 ppm)</i>   | < 3 mV/ day (8 hours)                                |
| <i>(Measured at constant temperature and with ISE and Reference Electrode continually immersed)</i>      |  |

## Interference:

All poly-crystalline membranes contain Silver Sulphide and thus will not give reliable readings if Ag or S ions are present in the solution. Cyanide has a Selectivity Coefficient (SC) of 1 (equally sensitive to CN and I) and hence will cause a significant positive error if it is present in concentrations greater than one one-hundredth of the iodide. Other minor interferences are from Br (SC~ 0.0004) and Cl (0.000001). The SC is the approximate apparent increase in the measured concentration caused by 1 unit of the interferent. Thus the likely effect of any interfering ion (% increase) can be calculated as follows:

$$((\text{expected concentration}) \times (\text{SC}) / (\text{expected Iodide concentration})) \times 100.$$

**For more information, see: [www.nico2000.net](http://www.nico2000.net).**